

MS of 23052-19-9

MASS SPECTROMETRY REPORT



Sample Description	Instrument	Waters ZQ2000
Analyzed date: 2024-04-10	Probe:	ESI
Analysis: PY	Nebulizer Gas Flow:	1.5L/min
Sample: 氧化 r-EC	CDL:	-20.0v
M.W.: 498.50	CDL Temp.:	250 °C
Lot. No.: 氧化 r-EC	Block Temp.:	200 °C
	Probe Bias:	+4.0kv
	Detector:	1.5kv
	T. Flow:	0.2ml/min
	B. Conc.:	50%H2O/50%ACN

*For laboratory or manufacturing use only, not for drug, food and household use.

* If you have any comments and suggestions, please contact us by sales@amadischem.com