


MS of 135861-77-7



Sample Description	Instrument	Waters ZQ2000	
Analyzed date	Probe	ESI	Probe Bias: +4.5kV
Analyst	Nebulizer Gas Flow	1.5L/min	Detector: 1.5kV
Sample	CDL	-20.0V	T. Flow: 0.2ml/min
M.W.	CDL Temp.	250°C	B. Conc.: 50%H2O 50%ACN
Lot No.	Block Temp.	200°C	

* If you have any comments and suggestions, please contact us by sales@amadischem.com